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11. Bjerring P. Photorejuvenation—An overview. *Med Laser Appl* 2004;19:186–195.
12. Ross EV, Barnette DJ, Glatter RD, Grevelink JM. Effects of overlap and pass number in CO₂ laser skin resurfacing: A study of residual thermal damage, cell death, and wound healing. *Lasers Surg Med* 1999;24(2):103–112.
13. Ross EV, McKinlay JR, Sajben FP, Miller CH, Barnette DJ, Meehan KJ, Chhieng NP, Deavers MJ, Zelickson BD. Use of a novel erbium laser in a Yucatan minipig: A study of residual thermal damage, ablation, and wound healing as a function of pulse duration. *Lasers Surg Med* 2002;30(2):93–100.
14. Majaron B, Srinivas SM, Huang H, Nelson JS. Deep coagulation of dermal collagen with repetitive Er:YAG laser irradiation. *Lasers Surg Med* 2000;26(2):215–222. 5
15. Majaron B, Verkruysse W, Kelly KM, Nelson JS. Er:YAG laser skin resurfacing using repetitive long-pulse exposure and cryogen spray cooling: II. Theoretical analysis. *Lasers Surg Med* 2001;28(2):131–137.
16. Kao B, Kelly KM, Majaron B, Nelson JS. Novel model for evaluation of epidermal preservation and dermal collagen remodeling following photorejuvenation of human skin. *Lasers Surg Med* 2003;32(2):115–119.
17. Drnovsek-Olup B, Beltram M, Pizem J. Novel model for evaluation of epidermal preservation and dermal collagen remodeling following photorejuvenation of human skin. *Lasers Surg Med* 2003;32(2):115–119.
18. Ross EV, Miller C, Meehan K, McKinlay J, Sajben P, Trafeli JP, Barnette DJ. One-pass CO₂ versus multiple-pass Er:YAG laser resurfacing in the treatment of rhytides: A comparison side-by-side study of pulsed CO₂ and Er:YAG lasers. *Dermatol Surg* 2001;27(8):709–715
19. Adrian RM. Pulsed carbon dioxide and long pulse 10-ms erbium-YAG laser resurfacing: A comparative clinical and histologic study. *J Cutan Laser Ther* 1999;1(4):197–202.
20. Newman JB, Lord JL, Ash K, McDaniel DH. Variable pulse erbium:YAG laser skin resurfacing of perioral rhytides and side-by-side comparison with carbon dioxide laser. *Lasers Surg Med* 2000;26(2):208–214.
21. Drosner M. Non-ablative wrinkle reduction by Er:YAG laser: Encouraging results (case report). *Lasers Surg Med* 2002; Suppl 14:66.
22. Kunzi-Rapp K, Cambier B, Drosner M, Dierickx C, Levy JL, Russ D, Steiner R. Non-ablative skin rejuvenation with Erbium:YAG laser pulses—investigation of structural changes in the skin. *Lasers Med Sci* 2003;18(Suppl 1):1–68.
23. Kunzi-Rapp K, Cambier B, Drosner M, Dierickx C, Levy JL, Larrouy JC, Russ D, Steiner R. Non-ablative versus ablative Er:YAG laser therapy in photoaged skin: Histopathological changes. *Lasers Surg Med* 2004; Suppl 16:21.
24. Drnovs'ek-Olup B, Beltram M, Piz'em J. Repetitive Er:YAG laser irradiation of human skin: A histological evaluation. *Lasers Surg Med* 2004;35:146–151.
25. Trelles MA, Allones I, Luna R. Facial rejuvenation with a nonablative 1320 nm Nd:YAG laser: A preliminary clinical and histologic evaluation. *Dermatol Surg* 2001;27(2):111–116.
26. Maneker GM, Wrone DA, WilliamRM, Moy RL. Treatment of facial rhytides with a non-ablative laser; a clinical and histologic study. *Dermatol Surg* 1999;25:440–444.
27. Stuzin JM, Baker TJ, Baker TM, Kligman AM. Histologic effects of the high-energy pulsed CO₂ laser on photoaged facial skin. *Plast Reconstr Surg* 1997;99(7):2036–
28. Kligman AM, Zheng P, Lavker RM. The anatomy and pathogenesis of wrinkles. *Br J Dermatol* 1985;113(1):37–42.
29. Fisher GJ, Kang S, Varani J, Bata-Csorgo Z, Wan Y, Datta S, Voorhees JJ. Mechanisms of photoaging and chronological skin aging. *Arch Dermatol* 2002;138(11):1462–1670.
30. El-Domyati M, Attia S, Saleh F, Brown D, Birk DE, Gasparro F, Ahmad H, Uitto J. Intrinsic aging vs. photoaging: A comparative histopathological, immunohistochemical, and ultrastructural study of skin. *Exp Dermatol* 2002;11(5):